


So easy to use... It's automatic!

- Auto-access: automatic aperture of the slot device via motion detector
- Auto-detect: **PSPIX** automatically recognizes the size of the inserted plate
- Auto-scan: automatic readout after insertion of the Imaging Plate (avg 5 sec.)
- Auto-erase: the Imaging Plate is automatically erased after scanning
- Auto-eject: the Imaging Plate is automatically ejected in its collector for immediate reuse
- Auto-optimization: **PSPIX** adjusts the brightness and contrast of the image automatically - Image quality is always optimized
- Auto-sleep: **PSPIX** automatically switches to the energy-saving mode when it is not used for a certain period of time

pspix

All the advantages of conventional film; flexibility and thinness... yet digital and reusable!



Available in sizes 0 through 3



100 %
ACTIVE AREA



Onboard display

ONE

You shoot

TWO

You scan




Auto detect..



Auto scan..



Auto eject!



sopro
Imaging



FREE!

* With purchase of PSPIX