

QUICKRAY[®]
Trusted Since 1998



The
HASSLE FREE
Intraoral Sensor
Low Cost, Never Any Hidden Fees

The QuickRay Digital X-Ray Sensor series gives you excellent image clarity at the highest resolution with minimum radiation exposure; and connects directly to your computer's USB port without requiring a special interface box.



TWAIN compatible, QuickRay works in all leading imaging software and can be used in conjunction with existing sensors.

COMPATIBLE WITH





Your Practice Deserves To Purchase The *Right* Sensor!

Let Our Team of Experts Help You Make the Best Choice

SENSOR TECHNICAL DATA

SENSOR TECH	CMOS
RESOLUTION	Line Pairs >20 lp/mm
ACTIVE AREA	Size 1: 19.95mm x 30.02mm Size 2: 25.99mm x 36.00mm
SCINTILLATOR	CSI
PIXEL SIZE	19 x19 Microns
COMPATIBILITY	Windows 32/64 Bit & MAC OS
BODY THICKNESS	5.3 MM



QUICKRAY SENSOR



- Proven Durable Engineering
- Hermetically Sealed Housing
- Shock Absorbing Material
- Ergonomic Design
- Thin with Rounded Corners for Patient Comfort

CORD



- 2M Cable Length
- Kevlar Reinforced Virtually Unbreakable
- Replaceable
- Direct High Speed 2.0 USB Connectivity



ACTUAL IMAGE

SENSOR KIT INCLUDES



Protective Sheaths



Rinn Positioning Kit



Wall Mount Holder



Warranty

2 Year Warranty
7 Year Replacement Plan



Expedited Rental Service

Available for NO Downtime!

CONTACT US TODAY!



Outstanding Dental Equipment Since 1989

www.VideoDental.com | sales@videodental.com | 800.323.2690

100 East Granada Boulevard Suite 219 Ormond Beach, FL 32176